MESTRA

No disponible para:
(C.E.E.)

CENTRÍFUGA G3 R-040000

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*Not available for:
(C.E.E.)*

CENTRIFUGE G3 R-040000

As with other models, MESTRA has in this case based its design on principles of strict safety. Features include an external fastening lever and firing outside the circle generated by the arm.

The arm of the MESTRA G3 centrifuge is articulated. This facilitates the casting of highly-melted materials by adding rotary speed to initial speed.


The counterweight can be displaced gradually according to the weight and size of the cylinder.

Being height-adjustable, the cylinder holder can accommodate any MESTRA cylinder within its crucible mouth.

The built-in crucible is 99% silicon oxide. This means more castings and higher purity.

The body of MESTRA G3 centrifuge is specially strong. Made from special alumina, it is modern and harmonious in appearance.

USE

1. Load clockwise, turning two and half to three times. Then secure by moving handle upward.
2. Counterbalance according to cylinder size.
3. Place the articulated arm in position. 
4. Place the preheated crucible on its support.
5. Load metal into the crucible and melt it. Apply heat to keep it melted for several minutes.
6. Once it has been preheated, place the cylinder on its support. Ensure that its orifice is as that of the crucible and that they are well fitted.
7. Start the centrifuge by pushing the handle down.

DIMENSIONS	
Diameter of turning circle	540 mm
Max.height	260 mm
Weight	13 kg

IMPORTANT

- Ensure that the screw fixing the shaft to the arm is tight.
- Ensure that the arm is properly counterweighted (see use).
- Retract the counterweight and the cylinder/crucible assembly, turning the two articulations towards the centre of the apparatus so that the whole assembly finally has the shape shown in the figure:

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