Hornos de PRECALENTAMIENTO

HP-100 R-080148

HP-50 R-080144

HP-25 R-080142

HP-100 R-080148

HP-50 R-080144

HP-25 R-080142

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The family of HP furnaces is made up by three models that only differ in their muffle capacity. Depending on the needs of each laboratory, it is possible to choose the most suitable furnace in order to optimize the energy expenditure in accordance to the casting rings quantity that is usually preheated. The range begins with the HP-25 (2,5 litres muffle capacity), an adequate option for reduce production laboratories. It continues with the HP-50 (5 litres muffle capacity), adequated for medium level production. The HP-100, with 10 litres muffle volume represents the highest capacity model, recommended for high production laboratories.

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
HP family muffle chassis is entirely made of stainless steel, heat degradation resistant material and easy to clean. It also has an excellent finish.

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The muffle is made in ceramic material with heaters wrapped along the upper, lower and on both sides. This way, there is a uniform interior distribution of temperatures. The exterior sides of the muffle are thermally insulated by zirconium silicate, probably the best thermal insulator used in industrial processes.

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Gas flow inside the muffle has been carefully studied. There is an air inlet at the door and a ventilation pipe at the back, both with calibrated sections. This way air inside the muffle is renovated without reducing thermal efficiency of the oven nor producing temperature differences.

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
A practical socket allows the synchronised operation of the oven together with an extraction unit. Automatic disconnection can also be programmed between 75° C and 850° C. When temperature is over 850° C, the control automatically disconnects the power supply.

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Magnetic closing system is easy and simple. A safety sensor detects door opening, and disconnects heaters while it is opened.

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☒ TECLA DE PROGRAMAS (20 Prog.)

ERGONOMICS

Twin display screen: one for temperature, the other for time. Graphic showing temperature increase and holding, with indication of what process is ongoing.

EASY TO HANDLE

- ▶ **SEARCH KEY:** press to enter programming mode. The setting to be changed on the graph flashes.
- ▲ **UP KEY:** press to increase the time or temperature setting flashing on the graphic.
- ▼ **DOWN KEY:** press to decrease the time or temperature setting flashing on the graphic.
- ☒ **ESCAPE KEY:** press to exit programming mode and memorise settings.
- ☒ **STOP/START KEY:** hold down for 5 seconds to turn the oven on and off.

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
With this electronic control run by a microprocessor temperature can be defined accurately in the three heating ramps and in the three plateaux that can be programmed in a cycle, and also the required times in order to reach these temperatures. As the most relevant characteristics of the control it can be pointed out the following.

- Continuous monitorization of temperatures and times through its two displays.
- Precise control of temperatures, with grade by grade signal.
- Heating resistance control by triac with digital switch on when going over "0".
- PID algorithm Software for a best control of thermal inertia.
- Memory of up to 20 programmes.
- Synoptic panel indicating the stage in which we are every time.
- Thermocouple breakage indicator.
- Warning buzzing of the most significant steps of the cycle.
- It is possible to programme a time relay switch on. (99 h. 59')
- Overvoltage or voltage variation protection.

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The Smoke Extractor R-080125 is an interesting attachment to the HP Preheating Furnaces. Installing this device, you succeed in getting rid of furnace gases generated in the muffle in a continuous and smooth way. It is completely made of stainless steel and its installation is easy. If it is connected to the MESTRA HP furnaces socket, the extractor ranges from 75° to 850°c as programmed by the user. If it is connected to another socket it runs continuously.

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