
Ref. 080620

Nº	Ref.	Concepto	Description
1	080620-10	Junta tapa	Lid O-ring
2	000334	Racor recto	Straight fitting
3	080620-04	Válvula entrada	Inlet valve
4	080620-03	Regulador presión	Pressure regulator
5	080620-13	Válvula seguridad	Safe valve
6	080620-12	Acoplador	Connector
7	080250-14	Manómetro	Pressure gauge
8	080620-05	Válvula salida	Outlet valve
9	080620-15	Tapa con válvula de apertura (sin junta)	Lid with open valve (without O-ring)

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Ignacio Mestraitua  
Director Gerente  
Sondika, 15.02.2021

Ref. 080620 is a tank specially designed for the storage and supply of distilled water under pressure in the dental sector.

## Before filling the tank

Before filling the tank, it must be properly cleaned and sterilized. Do not use detergents containing hydrochloric acid, hypochlorous acid, etc.

## Installation

For correct use of the appliance, we recommend that you adhere to the following instructions:

- Carefully unpack the appliance.
- Place the tank away from heat sources.
- Check that it has not suffered any damage during transport.
- Check that the tank is not scratched or dented.
- Close the input and output valves (B and E).
- Connect pressurized air to inlet connector (B).
- Connect the power tube for your equipment to the outlet connector (E).
- Open the lid and fill the tank with water.
- Close the lid.
- Open the pressurized air inlet cock (A).
- Use the regulator (D) to adjust the pressure (below 5 bar).
- Open the water outlet tap (E).

## Filling

- To open the tank, first close the inlet and outlet valves (B and E).
- Then, pull the ring of the opening valve (A) to depressurize the tank.
- Now you can open the lid and fill the tank.
- Close the lid.
- Open the inlet and outlet valves (B and E).

## Maintenance and cleaning

To ensure correct, reliable, and enduring operation, we recommend:

- Use the tank with care.
- Avoid falls, bumps and inappropriate movements of the tank; they could cause it to malfunction.
- Do not use if the tank is damaged, scratched or dented.
- Do not use if the cover is dented or bent.
- The cover and connections must always be kept clean and free of any contaminating agent.

- The life of the lid's O-ring depends on the use of the tank, the storage location and the environmental conditions. It is recommended that this component be replaced at least once a year.
- The tank should be located away from heat sources.

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
Meets the requirements set out in the directive:

2014/68/CE (pressure equipment)

Ignacio Mestraitua  
Managing Director  
Sondika, 2021-02-15